

Feasibility of Non-Invasive Monitoring of Respiratory, ECG and Activity Parameters in Pair-Housed Freely-Moving Adult Male Rats using a Jacketed External Telemetry System (DECRO®); A Pharmacological Validation with Baclofen.



S Baudet¹, R Fares², J Vial¹, M Ravaz¹, V Roger¹, C Bory¹, T Flénet²
¹ Charles River Laboratories France Safety Assessment SAS, Lyon, France
² Etisense SAS, Lyon, France

1 INTRODUCTION

Simultaneous assessment of multiple vital functions, along with collection of standard toxicology endpoints in non-rodents using non-invasive telemetry jackets, is increasingly conducted in regulatory non-clinical studies. Equivalent solutions are lacking in rats, despite being the most common rodent species used in toxicology studies.

A Bluetooth-based telemetry jacket (DECRO solution) was developed to non-invasively record respiratory function, ECG and activity in rats.

A feasibility study was conducted in group-housed rats in environmental/husbandry conditions similar to those employed in 4-week regulatory studies.

2 MATERIAL AND METHODS


ANIMALS AND TEST ARTICLES

Animals
 Pair-housed, freely-moving male Wistar rats
 Age: 9 weeks on first dosing occasion


Test materials
 Vehicle control (n=4)
 R(+)-Baclofen (15 mg/kg; n=4)

Gavage dosing

EQUIPMENT



EXPERIMENTAL DESIGN



PARAMETERS AND DATA ANALYSIS

- Proprietary software DECRO 1.0.0
- 10-minute averages of good quality, continuously or intermittently acquired telemetry signals (-2h to +24h)
- Recorded parameters: heart rate, respiratory rate, inspiratory/expiratory times, tidal/minute volumes, activity level (10-minute averages)

STATISTICS

- Two-way (time and dose) repeated measures ANOVA
- If significant time-by-dose interaction, vehicle vs baclofen comparisons at each time-point
- If no interaction, global dose effect (vehicle vs baclofen) averaged over all time-points.

4 DISCUSSION


- DECRO® external telemetry jacket was very well tolerated when worn for 24 hours on several occasions by pair-housed, growing rats
- The quality of the respiratory, ECG and accelerometry signals permitted to robustly analyze the effects of baclofen on these vital functions
- Baclofen (15 mg/kg, p.o.) exerted the expected depressant effects on respiration (lower respiratory rate/higher inspiratory time) and on the level of activity. Baclofen effects generally resumed by 3.5 to 6 hours postdose
- To our knowledge, this study revealed a novel finding: **the sustained and reproducible baclofen-induced increase in heart rate in freely-moving, conscious rats.**
- The tachycardic effects of baclofen occurred while the activity level was very low. This finding supports:
 - An intrinsic pharmacological effect of baclofen on cardiac function
 - The interest of assessing the animal's activity when testing new compounds

5 CONCLUSION

This study demonstrates the feasibility of non-invasive and longitudinal assessment of vital physiological functions using the DECRO telemetry jacketed solution in socially-housed, freely-moving rats, in both safety pharmacology and toxicology experimental conditions.

3 RESULTS / KEY OUTCOMES

HEALTH STATUS AND WELL-BEING




JACKET TOLERABILITY


- No general and local clinical signs with repeated/intermittent jacket wearing
- Normal rat growth with repeated/intermittent jacket wearing and baclofen administration
- Expression of normal behavioral repertoire and social interactions

DATA ACQUISITION AND QUALITY

- Social housing did not preclude to record a high load of good to excellent quality signal (> 95%)
- Signal filtering was done mostly based on exclusion of non-physiological values
- Markers for respiratory/ECG cycle identification were properly placed with the automated algorithm



EFFECTS OF BACLOFEN EVALUATED WITH DECRO JACKET



Week 1: acute assessment ("ICH S7A")
 Wks 2 and 4: longitudinal assessment ("ICH M3(R2)")

	Heart Rate	Respiratory Rate	Inspiratory Time	Expiratory Time	Tidal Volume	Minute Volume	Activity Level	Statistical Significance
Week 1	↗↗	↘↘	↘↘	↗	↔↔	↘	↘↘	1/1.5h to 4.5/6h
Week 2	↗	↘	↘	↔↔	↗	↔↔	↘	1/2h to 3.5/6h
Week 4	↗	↘	↘	↔↔	↔↔	↔↔	↔↔	1h to 6h